



# SLC30A10 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA747695LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6XR72
<b>Immunogen</b>	Recombinant Human Zinc transporter 10 protein (58-278AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Zinc transporter 10 (ZnT-10) (Manganese transporter SLC30A10) (Solute carrier family 30 member 10), SLC30A10, ZNT10 ZNT8
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC30A10
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.