



SLC30A10 Antibody, Biotin conjugated

Product Code	CSB-PA747695LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XR72
Immunogen	Recombinant Human Zinc transporter 10 protein (58-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 10 (ZnT-10) (Manganese transporter SLC30A10) (Solute carrier family 30 member 10), SLC30A10, ZNT10 ZNT8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC30A10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.