



WDFY4 Antibody, Biotin conjugated

Product Code	CSB-PA747798LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZS81
Immunogen	Recombinant Human WD repeat- and FYVE domain-containing protein 4 protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat- and FYVE domain-containing protein 4, WDFY4, C10orf64 KIAA1607
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDFY4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.