



FAAP20 Antibody, Biotin conjugated

Product Code	CSB-PA750774HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NZ36
Immunogen	Recombinant Human Fanconi anemia core complex-associated protein 20 protein (14-179AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fanconi anemia core complex-associated protein 20 (FANCA-associated protein of 20 kDa) (Fanconi anemia-associated protein of 20 kDa), FAAP20, C1orf86
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAAP20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.