



THADA Antibody, Biotin conjugated

Product Code	CSB-PA751092LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6YHU6
Immunogen	Recombinant Human Thyroid adenoma-associated protein (432-651AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroid adenoma-associated protein (Gene inducing thyroid adenomas protein), THADA, GITA KIAA1767
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	THADA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.