



ZNF365 Antibody, Biotin conjugated

Product Code	CSB-PA754628LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70YC5
Immunogen	Recombinant Human Protein ZNF365 protein (241-401AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ZNF365 (DISC1-binding zinc-finger protein) (Protein su48), ZNF365, DBZ KIAA0844
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF365
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.