



FGD5 Antibody, Biotin conjugated

Product Code	CSB-PA757939LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNL6
Immunogen	Recombinant Human FYVE, RhoGEF and PH domain-containing protein 5 protein (269-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FYVE, RhoGEF and PH domain-containing protein 5 (Zinc finger FYVE domain-containing protein 23), FGD5, ZFYVE23
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FGD5