




PRMT9 Antibody

Product Code	CSB-PA764718LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P2P2
Immunogen	Recombinant Human Putative protein arginine N-methyltransferase 9 protein (181-453AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein arginine N-methyltransferase 9 (Protein arginine N-methyltransferase 10) (EC 2.1.1.320), PRMT9, PRMT10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRMT9

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA764718LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.