




DTX4 Antibody

Product Code	CSB-PA764784LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2E6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase DTX4 protein (167-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase DTX4 (EC 2.3.2.27) (Protein deltex-4) (Deltex4) (RING finger protein 155) (RING-type E3 ubiquitin transferase DTX4), DTX4, KIAA0937 RNF155
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DTX4

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA764784LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.