



DTX4 Antibody, Biotin conjugated

Product Code	CSB-PA764784LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2E6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase DTX4 protein (167-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase DTX4 (EC 2.3.2.27) (Protein deltex-4) (Deltex4) (RING finger protein 155) (RING-type E3 ubiquitin transferase DTX4), DTX4, KIAA0937 RNF155
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DTX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.