



KANSL1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA770344LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7Z3B3 |
| Immunogen | Recombinant Human KAT8 regulatory NSL complex subunit 1 protein (815-1086AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | KAT8 regulatory NSL complex subunit 1 (MLL1/MLL complex subunit KANSL1) (MSL1 homolog 1) (hMSL1v1) (NSL complex protein NSL1) (Non-specific lethal 1 homolog), KANSL1, CENP-36 KIAA1267 MSL1V1 NSL1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | KANSL1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |