



SBF2 Antibody, HRP conjugated

Product Code	CSB-PA773047LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WG5
Immunogen	Recombinant Human Myotubularin-related protein 13 protein (1585-1813AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 13 (SET-binding factor 2), SBF2, CMT4B2 KIAA1766 MTMR13
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SBF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.