



## SBF2 Antibody, Biotin conjugated

Product Code	CSB-PA773047LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WG5
Immunogen	Recombinant Human Myotubularin-related protein 13 protein (1585-1813AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 13 (SET-binding factor 2), SBF2, CMT4B2 KIAA1766 MTMR13
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SBF2