



MILR1 Antibody, Biotin conjugated

Product Code	CSB-PA801819LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6M3
Immunogen	Recombinant Human Allergin-1 protein (249-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Allergin-1 (Allergy inhibitory receptor 1) (Mast cell antigen 32) (MCA-32) (Mast cell immunoglobulin-like receptor 1), MILR1, C17orf60 MCA32
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MILR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.