



MORC1 Antibody, Biotin conjugated

Product Code	CSB-PA803127LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VD1
Immunogen	Recombinant Human MORC family CW-type zinc finger protein 1 protein (683-984AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORC family CW-type zinc finger protein 1 (Cancer/testis antigen 33) (CT33), MORC1, MORC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MORC1