





PLIN5 Antibody

Product Code	CSB-PA804867LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00G26
Immunogen	Recombinant Human Perilipin-5 protein (305-453AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perilipin-5 (Lipid storage droplet protein 5), PLIN5, LSDP5 OXPAT PAT-1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLIN5

Image



Western blot
All lanes: PLIN5 antibody at 2 µg/ml
Lane 1: Mouse kidney tissue
Lane 2: Mouse brain tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 51 kDa
Observed band size: 51, 36 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA804867LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

