



PLIN5 Antibody, HRP conjugated

Product Code	CSB-PA804867LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00G26
Immunogen	Recombinant Human Perilipin-5 protein (305-453AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perilipin-5 (Lipid storage droplet protein 5), PLIN5, LSDP5 OXPAT PAT-1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLIN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.