




ALDH16A1 Antibody

Product Code	CSB-PA811634LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZ83
Immunogen	Recombinant Human Aldehyde dehydrogenase family 16 member A1 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aldehyde dehydrogenase family 16 member A1, ALDH16A1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ALDH16A1

Image



Immunofluorescent analysis of HeLa cells using CSB-PA811634LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.