



TP53INP2 Antibody, HRP conjugated

Product Code	CSB-PA812888LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXH6
Immunogen	Recombinant Human Tumor protein p53-inducible nuclear protein 2 protein (100-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible nuclear protein 2 (Diabetes and obesity-regulated gene) (p53-inducible protein U) (PIG-U), TP53INP2, C20orf110 DOR PINH
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53INP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.