



FMNL3 Antibody, FITC conjugated

Product Code	CSB-PA815554LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVF7
Immunogen	Recombinant Human Formin-like protein 3 protein (2-212AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Formin-like protein 3 (Formin homology 2 domain-containing protein 3) (WW domain-binding protein 3) (WBP-3), FMNL3, FHOD3 FRL2 KIAA2014 WBP3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FMNL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.