



P3H2 Antibody, Biotin conjugated

Product Code	CSB-PA815555HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVL5
Immunogen	Recombinant Human Prolyl 3-hydroxylase 2 protein (301-528AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolyl 3-hydroxylase 2 (EC 1.14.11.7) (Lepreca-like protein 1) (Myxoid liposarcoma-associated protein 4), P3H2, LEPREL1 MLAT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	P3H2