



DOCK3 Antibody, Biotin conjugated

Product Code	CSB-PA818270LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZD9
Immunogen	Recombinant Human Deducator of cytokinesis protein 3 protein (1630-1934AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deducator of cytokinesis protein 3 (Modifier of cell adhesion) (Presenilin-binding protein) (PBP), DOCK3, KIAA0299 MOCA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DOCK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.