



DMXL2 Antibody, FITC conjugated

Product Code	CSB-PA819469LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDJ6
Immunogen	Recombinant Human DmX-like protein 2 protein (2471-2657AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DmX-like protein 2 (Rabconnectin-3), DMXL2, KIAA0856
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DMXL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.