



RTF1 Antibody, Biotin conjugated

Product Code	CSB-PA821675HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92541
Immunogen	Recombinant Human RNA polymerase-associated protein RTF1 homolog protein (620-670AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA polymerase-associated protein RTF1 homolog, RTF1, KIAA0252
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RTF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.