



# SIRT1 Antibody, HRP conjugated

|                            |                                                                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA822202HB01HU                                                                                                                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                       |
| <b>Uniprot No.</b>         | Q96EB6                                                                                                                                                                                              |
| <b>Immunogen</b>           | Recombinant Human NAD-dependent protein deacetylase sirtuin-1 protein (2-747AA)                                                                                                                     |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                              |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                               |
| <b>Form</b>                | Liquid                                                                                                                                                                                              |
| <b>Conjugate</b>           | HRP                                                                                                                                                                                                 |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                                                                    |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                            |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                          |
| <b>Alias</b>               | NAD-dependent protein deacetylase sirtuin-1 (hSIRT1) (EC 3.5.1.-) (Regulatory protein SIR2 homolog 1) (SIR2-like protein 1) (hSIR2) [Cleaved into: SirtT1 75 kDa fragment (75SirT1)], SIRT1, SIR2L1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                |
| <b>Research Area</b>       | Cell Biology                                                                                                                                                                                        |
| <b>Target Names</b>        | SIRT1                                                                                                                                                                                               |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                         |