



SIRT1 Antibody, FITC conjugated

Product Code	CSB-PA822202HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EB6
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-1 protein (2-747AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-1 (hSIRT1) (EC 3.5.1.-) (Regulatory protein SIR2 homolog 1) (SIR2-like protein 1) (hSIR2) [Cleaved into: SirtT1 75 kDa fragment (75SirT1)], SIRT1, SIR2L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SIRT1