



PDLIM5 Antibody, Biotin conjugated

Product Code	CSB-PA822223LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HC4
Immunogen	Recombinant Human PDZ and LIM domain protein 5 protein (268-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 5 (Enigma homolog) (Enigma-like PDZ and LIM domains protein), PDLIM5, ENH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PDLIM5