



PDLIM5 Antibody, Biotin conjugated

Product Code	CSB-PA822223LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HC4
Immunogen	Recombinant Human PDZ and LIM domain protein 5 protein (268-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 5 (Enigma homolog) (Enigma-like PDZ and LIM domains protein), PDLIM5, ENH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PDLIM5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.