



MEGF10 Antibody, Biotin conjugated

Product Code	CSB-PA822242LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KG7
Immunogen	Recombinant Human Multiple epidermal growth factor-like domains protein 10 protein (123-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multiple epidermal growth factor-like domains protein 10 (Multiple EGF-like domains protein 10), MEGF10, KIAA1780
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MEGF10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.