



PRIMPOL Antibody, FITC conjugated

Product Code	CSB-PA822255LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LW4
Immunogen	Recombinant Human DNA-directed primase/polymerase protein (305-537AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed primase/polymerase protein (hPrimpol1) (EC 2.7.7.-) (Coiled-coil domain-containing protein 111), PRIMPOL, CCDC111
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRIMPOL