



PRIMPOL Antibody, Biotin conjugated

Product Code	CSB-PA822255LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LW4
Immunogen	Recombinant Human DNA-directed primase/polymerase protein (305-537AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed primase/polymerase protein (hPrimpol1) (EC 2.7.7.-) (Coiled-coil domain-containing protein 111), PRIMPOL, CCDC111
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRIMPOL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.