



TSNARE1 Antibody, HRP conjugated

Product Code	CSB-PA822269LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NA8
Immunogen	Recombinant Human t-SNARE domain-containing protein 1 protein (6-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	t-SNARE domain-containing protein 1, TSNARE1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSNARE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.