




Nlrp3 Antibody

Product Code	CSB-PA823181LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R4B8
Immunogen	Recombinant Mouse NACHT, LRR and PYD domains-containing protein 3 protein (1-153AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 3 (Cold autoinflammatory syndrome 1 protein homolog) (Cryopyrin) (Mast cell maturation-associated-inducible protein 1) (PYRIN-containing APAF1-like protein 1), Nlrp3, Cias1 Mmig1 Nalp3 Pypaf1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Nlrp3

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA823181LA01MO at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.