



# Nlrp3 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA823181LD01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8R4B8
<b>Immunogen</b>	Recombinant Mouse NACHT, LRR and PYD domains-containing protein 3 protein (1-153AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	NACHT, LRR and PYD domains-containing protein 3 (Cold autoinflammatory syndrome 1 protein homolog) (Cryopyrin) (Mast cell maturation-associated-inducible protein 1) (PYRIN-containing APAF1-like protein 1), Nlrp3, Cias1 Mmig1 Nalp3 Pypaf1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Immunology
<b>Target Names</b>	Nlrp3