



Nlrp3 Antibody, Biotin conjugated

Product Code	CSB-PA823181LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R4B8
Immunogen	Recombinant Mouse NACHT, LRR and PYD domains-containing protein 3 protein (1-153AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 3 (Cold autoinflammatory syndrome 1 protein homolog) (Cryopyrin) (Mast cell maturation-associated-inducible protein 1) (PYRIN-containing APAF1-like protein 1), Nlrp3, Cias1 Mmig1 Nalp3 Pypaf1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Nlrp3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.