



P4HB Antibody

Product Code	CSB-PA823188YA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R4U2
Immunogen	Recombinant Cricetulus griseus Protein disulfide-isomerase protein (20-509AA)
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Prolyl 4-hydroxylase subunit beta) (p58), P4HB, PDIA1
Immunogen Species	Cricetulus griseus
Research Area	Others
Target Names	P4HB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.