



RNF125 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA836205LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96EQ8 |
| Immunogen | Recombinant Human E3 ubiquitin-protein ligase RNF125 protein (143-231AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | E3 ubiquitin-protein ligase RNF125 (EC 2.3.2.27) (RING finger protein 125) (T-cell RING activation protein 1) (TRAC-1), RNF125 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | RNF125 |