



# RNF125 Antibody, Biotin conjugated

Product Code	CSB-PA836205LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EQ8
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF125 protein (143-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF125 (EC 2.3.2.27) (RING finger protein 125) (T-cell RING activation protein 1) (TRAC-1), RNF125
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF125