



PLD6 Antibody, HRP conjugated

Product Code	CSB-PA836649LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2A8
Immunogen	Recombinant Human Mitochondrial cardiolipin hydrolase protein (1-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial cardiolipin hydrolase (EC 3.1.-.-) (Choline phosphatase 6) (Mitochondrial phospholipase) (MitoPLD) (Phosphatidylcholine-hydrolyzing phospholipase D6) (Phospholipase D6) (PLD 6) (Protein zucchini homolog), PLD6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.