




C6orf223 Antibody

Product Code	CSB-PA836654LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N319
Immunogen	Recombinant Human Uncharacterized protein C6orf223 protein (9-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C6orf223, C6orf223
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C6orf223

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA836654LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.