



SPC24 Antibody, FITC conjugated

Product Code	CSB-PA836712LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBT2
Immunogen	Recombinant Human Kinetochores protein Spc24 protein (21-97AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetochores protein Spc24 (hSpc24), SPC24, SPBC24
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPC24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.