



SLC38A9 Antibody, HRP conjugated

Product Code	CSB-PA836713LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBW4
Immunogen	Recombinant Human Sodium-coupled neutral amino acid transporter 9 protein (1-119AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-coupled neutral amino acid transporter 9 (Solute carrier family 38 member 9) (Up-regulated in lung cancer 11), SLC38A9, URLC11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC38A9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.