



WNK4 Antibody, Biotin conjugated

Product Code	CSB-PA842690LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96J92
Immunogen	Recombinant Human Serine/threonine-protein kinase WNK4 protein (828-1106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase WNK4 (EC 2.7.11.1) (Protein kinase lysine-deficient 4) (Protein kinase with no lysine 4), WNK4, PRKWINK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WNK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.