



EXPH5 Antibody, Biotin conjugated

Product Code	CSB-PA843342LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEV8
Immunogen	Recombinant Human Exophilin-5 protein (1235-1490AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exophilin-5 (Synaptotagmin-like protein homolog lacking C2 domains b) (Slac2-b) (Slp homolog lacking C2 domains b), EXPH5, KIAA0624 SLAC2B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXPH5