



THRSP Antibody, Biotin conjugated

Product Code	CSB-PA846095LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92748
Immunogen	Recombinant Human Thyroid hormone-inducible hepatic protein (41-144AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroid hormone-inducible hepatic protein (Spot 14 protein) (S14) (SPOT14), THRSP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	THRSP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.