



FBXW5 Antibody, FITC conjugated

Product Code	CSB-PA846577LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969U6
Immunogen	Recombinant Human F-box/WD repeat-containing protein 5 protein (265-494AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 5 (F-box and WD-40 domain-containing protein 5), FBXW5, FBW5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.