



GALNT11 Antibody, Biotin conjugated

Product Code	CSB-PA847687LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCW6
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 11 protein (420-606AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polypeptide N-acetylgalactosaminyltransferase 11 (EC 2.4.1.41) (Polypeptide GalNAc transferase 11) (GalNAc-T11) (pp-GaNTase 11) (Protein-UDP acetylgalactosaminyltransferase 11) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 11), GALNT11
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GALNT11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.