




IL4I1 Antibody

Product Code	CSB-PA850426LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RQ9
Immunogen	Recombinant Human L-amino-acid oxidase protein (139-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IL4I1

Image



Western Blot
Positive WB detected in: HL-60 whole cell lysate(30μg)
All lanes: IL4I1 antibody at 1:1000
Secondary
Goat polyclonal to rabbit IgG at 1/40000 dilution
Predicted band size: 62,65 kDa
Observed band size: 70 kDa
Exposure time: 40s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.