



KCNK16 Antibody, Biotin conjugated

Product Code	CSB-PA850434LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96T55
Immunogen	Recombinant Human Potassium channel subfamily K member 16 protein (259-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 16 (2P domain potassium channel Talk-1) (TWIK-related alkaline pH-activated K(+) channel 1) (TALK-1), KCNK16, TALK1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KCNK16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.