



DNAAF3 Antibody, HRP conjugated

Product Code	CSB-PA850845LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9W5
Immunogen	Recombinant Human Dynein assembly factor 3, axonemal protein (1-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein assembly factor 3, axonemal, DNAAF3, C19orf51
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNAAF3